


MODEL 10 PART NUMBER MATRIX



Model Number:	10	Model 10		
Voltage:	1	120VAC, 60Hz	2	220VAC, 50/60Hz (CE Marked)
Plug Type:	B	Euro/Shuko*	I	Italy*
	C	Australia*	K	North America (Plug-in Ballast)**
	D	United Kingdom*	K	North America (In-Line Ballast)*
	E	Denmark*	L	Switzerland*
	F	France/Belgium*	R	Russia
	G	India*	S	Special (NEMA 6-20P)
	H	Israel*	X	No Plug (In-Line Style Only)
	Z	North America (NEMA 6-15P)		
	**Only plug type "K" is a plug-in style ballast, all other plug types are in-line ballast			
Adapter #:	#101	101	#122	122
	#102	102	#123	123
	#103	103	#125	#102 w/ 3" length
	#104	104	#127	#101 w/ 3" length
	#105	105	#128	#103 w/ 3" length
	#107	107	#206	206
	#108	108	#208	208
	#109	109	#211	211
	#110	110	#212	212
	#111	111	#213	213
	#118	118	#214	214
	#120	120	#216	216
			#217	217
				#218
				218
				#219
				219
				#221
				221
				#224
				224
				#229
				206 w/3" L
				#230
				208 w/3" L
				#231
				231
				#232
				232
				#233
				213 w/3" L
				#234
				234
				#236
				236
				#237
				237
				#238
				238
				#239
				239
				#240
				240
	#100	SM NO Adapter/Grommet	#115	SM Grommet
	#200	LG NO Adapter/Grommet	#215	LG Grommet
	* Adapter is threaded on both ends			
	Note: all adapters are available w/ ESD coating			
Lamp Type:	32	Warm White, 3200K	BL	Blue, 450nm
	42	Cool White, 4200K	GN	Green, 544nm
	51	White, 5100K	YL	Yellow, 585nm
	55	White, 5500K	RD	Red, 658nm
	BB	Black Light Blue, 365nm		
	Note: All White lamps are available w/ ESD coating			
Modifications:	t.b.d. (Consult KSE-Texas for available options or customization)			

ESD Material (Standard leave blank)
E1 for ESD



Microscope Adapter Cross Reference

#211	Aus Jena - Citival 2 Stereoscope OBS 330	#206	Leica MZ8	#110	Olympus SY3
#211	Clearview	#206	Leica MZ12.5	#211	Olympus SY4
#221	Denerica	#240	Leica S4E, S6E (male thrd, 1")	#103	Olympus SZ
#232	DVT	#206	Leica S4E	#236	Olympus SZ, (male/female thrd, 103 comparable)
#224	Fisher Stereomaster II	#118	Lecia SZ4	#206	Olympus SZX, SZH
#208	Fisher Sterezoom	#212	Meiji EMF, EMT 1, EMT 2	#208	Olympus VM
#108	Leica GZ4 (Male/female thread - 101 comparable)	#218	Meiji EMT 3	#110	Optem Zoom 65 Systems
#206	Leica GZ6	#219	Meiji EMX	#217	Parco, Meiji EMZ5
#231	Leica GZ6 (threaded version w/no objective)	#213	Meiji EMZ	#206	Prior
#104	Lecia (Bausch & Lomb) Academic Stereo	#212	Mitutoyo 377, Scope #12	#110	Science Scope
#101	Lecia (Bausch & Lomb) Stereo Zoom 1,2,3,4,5, STZ	#123	Mitutoyo TM510 2x Mag	#107	StockerYale Tool Analyzer
#118	Lecia (Bausch & Lomb) #6	#213	Mitutoyo Stereo	#109	Swift
#105	Lecia (Bausch & Lomb) #7	#102	Nikon SM5, SMZ2, SMZ6, SM6	#211	Unitron ZSB, ZST, FSB
#110	Lecia (American Optical) Cycloptic	#103	Nikon SMZ1, SMZ2B, SMZT2	#213	Unitron ZSM
#103	Lecia (American Optical) Stereo Star, male thrd 1"	#211	Nikon SMZ10	#206	Unitron (Motic)
#128	Lecia (American Optical) Stereo Star, male thrd 3"	#211	Nikon SMZ10	#206	Wild MSA, M3, M7, M8, MZ8, MZAPO
#110	Lecia Monozoom 7	#206	Nikon SMZ-U, 10A, 660, 645, 880	#206	Zeiss
		#237	Nikon SMZ660, SMZ645, male thrd, 1" long		
		#238	Nikon SMZ660, SMZ645, male thrd 3" long		
		#239	Nikon SMZ800, male thrd, 1" long		
		#211	Nissho or Bristol		
		#213	Nova Stereo		

